


# ECO Inject - Injection Needles

Single-use Injection Needle to Meet your Daily Routine



# ECO Inject - Injection Needles

Designed for Precise Injection, Safe Operation and Convenient Product Selection




103664




## Safety

Integrated safety spiral at the distal end prevents unintentional piercing through the sheath.




## Quality

Designed in four needle diameters and three needle lengths to offer what is needed for the daily routine.



## Time

Clear color codes allow a fast needle selection and secure a focused procedure.



## Fit


Three different working lengths enables an operation in the whole GI tract.

# ECO Inject - Injection Needles

**Designed for Precise Injection, Safe Operation and Convenient Product Selection**


## **ECO Inject**

Single-Use Injection Needle - 19G



## **ECO Inject**

Single-Use Injection Needle - 22G



## **ECO Inject**

Single-Use Injection Needle - 23G



## **ECO Inject**

Single-Use Injection Needle - 25G



# ECO Inject - Injection Needles

## Excellent Devices for a Wide Range of Specifications

### Ordering Information

Order No.	Article	Working length	Needle bevel	Needle length	Needle Ø	Quantity
E0420593	IN-O18-19-4-180	1800 mm	Regular	4 mm	19G	10
E0420594	IN-O18-19-5-180	1800 mm	Regular	5 mm	19G	10
E0420595	IN-O18-19-6-180	1800 mm	Regular	6 mm	19G	10
E0420596	IN-O18-22-4-180	1800 mm	Regular	4 mm	22G	10
E0420597	IN-O18-22-5-180	1800 mm	Regular	5 mm	22G	10
E0420598	IN-O18-22-6-180	1800 mm	Regular	6 mm	22G	10
E0420599	IN-O18-23-4-180	1800 mm	Regular	4 mm	23G	10
E0420600	IN-O18-23-5-180	1800 mm	Regular	5 mm	23G	10
E0420601	IN-O18-23-6-180	1800 mm	Regular	6 mm	23G	10
E0420602	IN-O18-25-4-180	1800 mm	Regular	4 mm	25G	10
E0420603	IN-O18-25-5-180	1800 mm	Regular	5 mm	25G	10
E0420604	IN-O18-25-6-180	1800 mm	Regular	6 mm	25G	10
E0420605	IN-O23-22-4-160	1600 mm	Regular	4 mm	22G	10
E0420606	IN-O23-22-5-160	1600 mm	Regular	5 mm	22G	10
E0420607	IN-O23-23-4-160	1600 mm	Regular	4 mm	23G	10
E0420608	IN-O23-23-5-160	1600 mm	Regular	5 mm	23G	10
E0420609	IN-O23-19-4-230	2300 mm	Regular	4 mm	19G	10
E0420610	IN-O23-19-5-230	2300 mm	Regular	5 mm	19G	10
E0420611	IN-O23-19-6-230	2300 mm	Regular	6 mm	19G	10
E0420612	IN-O23-22-4-230	2300 mm	Regular	4 mm	22G	10
E0420613	IN-O23-22-5-230	2300 mm	Regular	5 mm	22G	10
E0420614	IN-O23-22-6-230	2300 mm	Regular	6 mm	22G	10
E0420615	IN-O23-23-4-230	2300 mm	Regular	4 mm	23G	10
E0420616	IN-O23-23-5-230	2300 mm	Regular	5 mm	23G	10
E0420617	IN-O23-23-6-230	2300 mm	Regular	6 mm	23G	10
E0420618	IN-O23-25-4-230	2300 mm	Regular	4 mm	25G	10
E0420619	IN-O23-25-5-230	2300 mm	Regular	5 mm	25G	10
E0420620	IN-O23-25-6-230	2300 mm	Regular	6 mm	25G	10

 **[www.et-catalogue.olympus.eu](http://www.et-catalogue.olympus.eu)**

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications, accessories and service offerings may be required.

**OLYMPUS**

**OLYMPUS EUROPA SE & CO. KG**

Postbox 10 49 08, 20034 Hamburg, Germany

Wendenstrasse 20, 20097 Hamburg, Germany

Phone: +49 40 23773-0, Fax: +49 40 233765

[www.olympus-europa.com](http://www.olympus-europa.com)